



Microscope Wizards Ltd Risk Assessment & COSHH

Date:	Assessed by:	Position:	Date of last review:
24/04/2017	Alan J Monk	Senior Microscopy Specialist & Engineer	06/01/2025
<p>Task: Microscope servicing.</p> <p>A service is defined as a complete strip down of the optical components of the microscope which are then cleaned, reassembled (with lubrication as appropriate) and the optical and illumination path are realigned.</p> <p>Microscope Wizards Ltd is not responsible for the electrical safety testing of customer equipment and does not routinely provide Portable Appliance Testing (PAT) as part of the service. A visual inspection of the electrical components will be made by the engineer before the microscope is plugged in to the mains electricity supply (for the purposes of checking the bulb and alignment of the illumination pathway). If a fault is identified upon visual inspection, then the service will be terminated, and the customer informed and advised to take the defective equipment out of use immediately.</p> <p>Premises: Microscope Wizards Ltd workshop or customer premises (as defined by the job/customer).</p>			



Microscope Wizards Ltd Risk Assessment & COSHH

Risk Assessment					
Activity	Hazard	Who might be harmed and how	Existing measures to control risk	Risk rating	Result
Working within a laboratory environment – general hazards	Chemical	Engineer	When working at a customer's premises, Engineer should enquire where would be a safe place to perform servicing role to avoid being in close proximity with 'wet' laboratory work. Engineer must wear lab coat and nitrile gloves as minimum PPE when working in a laboratory setting.	Low	A
Movements and Postures	Musculoskeletal injury	Engineer, musculoskeletal pain/injury	Engineer to find a working position that is comfortable and to be aware of exercising good ergonomics at all times during the working day (avoiding stooping, over-reaching, twisting, fixed posture and ensuring there is adequate leg space when seated at a bench/desk). Engineer to avoid spending long periods, standing or sitting, alternate between postures as required to ensure comfort.	Low	A
Manual handling	Lifting, stretching, carrying	Engineer, musculoskeletal pain/injury	Engineer must not lift beyond personal capabilities and exercise good manual handling techniques at all times. Trolley to be used to assist when transporting heavy items to and from premises and vehicle.	Low	A
	Use of hand tools – screwdrivers and Allen keys	Engineer, musculoskeletal pain/injury	Tools should be correct for the application and of a good ergonomic design. Appropriate care should be taken to avoid repetitive movements and work should stop if discomfort is experienced.	Low	T
Connecting microscope to mains electricity	Electrical – risk of electrocution	Engineer	Engineer to perform a thorough visual inspection of the electrical components <u>before</u> plugging the microscope in to the mains electricity supply. If a fault is identified upon visual inspection or there are any signs of defect after connecting e.g., sparks then the service will be terminated and the customer advised to take the defective equipment out of use immediately.	Low	A


Microscope Wizards Ltd Risk Assessment & COSHH

Risk Assessment					
Activity	Hazard	Who might be harmed and how	Existing measures to control risk	Risk rating	Result
Application of cleaning solutions to microscope components	Chemical – risk to eyes and skin	Engineer	<p>Gloves to be worn to avoid risk of contact dermatitis associated with cleaning solutions.</p> <p>Cleaning solutions used are household grade only and comprise of foam cleaner and glass cleaner.</p> <p>Adequate eye protection to be worn if chemicals are to be sprayed.</p>	Low	T
Application of lubrication compounds to microscope components	Chemical	Engineer	<p>Gloves to be worn to avoid risk of contact dermatitis associated with lubricants.</p> <p>Lubricants used are standard mechanical greases and comprise of 3-in-1 oil, multi-purpose grease, silicone grease, graphite grease and WD40.</p> <p>Adequate eye protection to be worn if chemicals are to be sprayed.</p>	Low	T
Application of solvents to microscope components	<p>Chemical</p> <p>Flammables - solvents</p>	Engineer, others in the vicinity	<p>Gloves to be worn when handling solvents and work must be completed in a well-ventilated room.</p> <p>Due to the very small quantities involved it is not necessary to work within a fume hood.</p> <p>Solvents used: ethanol, acetone, IMS.</p> <p>Solvents are flammable – do not use near open flame or in close proximity to other potential sources of ignition. Small volumes significantly reduce this risk.</p>	Low	A



Microscope Wizards Ltd Risk Assessment & COSHH

COSHH								
Chemical	Hazard Symbol(s)	Signal Word	Quantity	Hazard Details	PPE & Other Safety Precautions	First Aid Measures	Risk Rating	Result
Ethanol, 96% laboratory grade	 	Danger	100ml	<ul style="list-style-type: none"> Highly flammable liquid & vapour Causes severe eye irritations 	<ul style="list-style-type: none"> Lab coat, gloves and eye protection must be worn Work within a well-ventilated area No open ignition sources 	<p>Skin contact:</p> <ul style="list-style-type: none"> Wash contaminated skin promptly with soap or mild detergent and water Remove clothing promptly, if soaked through, and wash as above <p>Eye contact:</p> <ul style="list-style-type: none"> Wash eyes promptly with plenty of water, while lifting the eye lids Continue to rinse for at least 15 minutes, and get medical attention <p>Ingestion:</p> <ul style="list-style-type: none"> Wash out mouth thoroughly, and give plenty of water to drink Seek medical attention <p>Inhalation:</p> <ul style="list-style-type: none"> Very small volumes in use Remove to fresh air 	Low	A

Microscope Wizards Ltd Risk Assessment & COSHH

COSHH								
Chemical	Hazard Symbol(s)	Signal Word	Quantity	Hazard Details	PPE & Other Safety Precautions	First Aid Measures	Risk Rating	Result
IMS, Industrial Methylated Spirit		Danger	100ml	<ul style="list-style-type: none"> Highly flammable liquid & vapour Causes severe eye irritations Harmful if swallowed May cause damage to organs 	<ul style="list-style-type: none"> Lab coat, gloves and eye protection must be worn Work within a well-ventilated area No open ignition sources 	<p>Skin contact:</p> <ul style="list-style-type: none"> Wash contaminated skin promptly with soap or mild detergent and water Remove clothing promptly, if soaked through, and wash as above <p>Eye contact:</p> <ul style="list-style-type: none"> Wash eyes promptly with plenty of water, while lifting the eye lids Continue to rinse for at least 15 minutes, and get medical attention <p>Ingestion:</p> <ul style="list-style-type: none"> Wash out mouth thoroughly, and give plenty of water to drink Seek medical attention <p>Inhalation:</p> <ul style="list-style-type: none"> Very small volumes in use Remove to fresh air 	Low	A

Microscope Wizards Ltd Risk Assessment & COSHH



COSHH								
Chemical	Hazard Symbol(s)	Signal Word	Quantity	Hazard Details	PPE & Other Safety Precautions	First Aid Measures	Risk Rating	Result
Acetone, analysis grade	 	Danger	100ml	<ul style="list-style-type: none"> Highly flammable liquid & vapour Causes severe eye irritations May cause drowsiness or dizziness Repeated exposure may cause skin dryness or cracking 	<ul style="list-style-type: none"> Lab coat, gloves and eye protection must be worn Work within a well-ventilated area No open ignition sources 	<p>Skin contact:</p> <ul style="list-style-type: none"> Wash contaminated skin promptly with soap or mild detergent and water Remove clothing promptly, if soaked through, and wash as above <p>Eye contact:</p> <ul style="list-style-type: none"> Wash eyes promptly with plenty of water, while lifting the eye lids Continue to rinse for at least 15 minutes, and get medical attention <p>Ingestion:</p> <ul style="list-style-type: none"> Wash out mouth thoroughly, and give plenty of water to drink Seek medical attention <p>Inhalation:</p> <ul style="list-style-type: none"> Very small volumes in use Remove to fresh air 	Low	A

Microscope Wizards Ltd Risk Assessment & COSHH

Action plan				
Ref No	Further action required	Action by whom	Action by when	Done
1				
2				
3				
4				
5				

Declaration by User

I confirm that I have read this Risk Assessment & COSHH and understand the hazards and risks involved. I confirm that I will follow all of the safety procedures as stated.

Name (printed): Alan J Monk	Signature: 	Date: 06/01/2025
Name (printed): Christine Monk	Signature: 	Date: 06/01/2025

Microscope Wizards Ltd Risk Assessment & COSHH

Risk Ratings:

LOW	unlikely that harm would arise under the controlled conditions listed, and even if exposure occurred, the injury would be relatively slight.
MEDIUM	more likely that harm might actually occur, and the outcome could be more serious e.g., some time off work, or a minor physical injury.
HIGH	an injury is likely to arise, and that injury might be serious (broken bones, trip to the hospital, loss of consciousness), or even a fatality.

Result:

T **Trivial risk.** The activity has been correctly identified as a hazard, but the risk is insignificant.

A **Adequately controlled, no further action necessary.** Suitable control measures are in place and the risk is low.

For T and A results, the assessment is complete.

N **Not adequately controlled, further action required.** The risk is medium or high, but it is reasonably practicable to reduce the risks further e.g., by purchasing mechanical aids for lifting. An action plan is required to state what needs to be carried out, by whom and within a suitable timeframe. Once complete, each line of the action plan must be signed off.